

Contents

4 Mount Type

6 Mount Type Options

8 Cart Type

10 Cart Type Options

12 Upholstery Colors

13 Space Requirement

14 Specifications







Flexible Assistant Table

The Assistant Table moves in a wide range to provide enough space for perfect assistance.

This table has

4 instrument holders.

1 High Volume Ejector

1 Saliva Ejector

3-way Syringe

1 Empty





Chair

- ① Hydraulic Chair
- ② Single Layer Cushion Backrest
- ③ Double Joint Headres
- 4 Left Armres
- ⑤ Single Layer Cushion Sea

Unit

- 6 Unit Body
- ② Assistant Table (4 Instrument holders
- ® Cuspidor/Spittoon Bowl (Glass)
- 9) Tissue Box Holder
- Water Temperature Indicator (Water heating system included
- ① Light Pol
- Dual Light ML 100A (Camera is not included

Delivery Table

- ① Delivery Table
- 🕦 Mouse Table
- (15) Secondary Tray
- **6** Control Pane
- 🛈 Instrument Holder (5 Holders
- (8) Foot Control
- ② 2nd Assi (Hands free suction)
- Micor NX Electric Motor Set (Instrument 3)
- ② Dmetec Scale

Stool

Eunctional Doctor Stock




















Chair

- ① Hydraulic Chair
- ② Single Layer Cushion Backrest
- ③ Double Joint Headres
- 4 Left Armres
- **(5)** Single Layer Cushion Seat

Unit


- 6 Unit Body
- ② Assistant Table (4 Instrument holders
- ® Cuspidor/Spittoon Bowl (Glass)
- 9) Tissue Box Holde
- Water Temperature Indicator(Water heating system included)
- 🛈 Light Pole
- Dual Light ML 100A (Camera is not included

Delivery Table

- ① Delivery Table
- 🕦 Mouse Table
- [®] Secondary Tray
- 16 Control Pane
- 🕏 Instrument Holder (5 Holders
- ® Foot Contro
- ② 2nd Assi (Hands free suction
- Micor NX Electric Motor Set (Instrument 3)
- ② Dmetec Scale


Stool

@ Functional Doctor Stoc











Upholstery Colors




For more information and AR experience, please visit imegagen.com.


Space Requirement



^{*} Minimum Space Requirement

Specifications







5	Unit	Body	_	
6	Dual Light ML100A (C	Camera is not included)		
7	Dbovine 5 Typ	e Side Table Set		
8	Assistant Table			
9	Water Temperature Indicator (Water heating system included)			
10	Cuspidor/Spittoon Bowl (Glass)			
11	Assistant Table (4 instrument holders): 1 High Vacuum Ejector, 1 Saliva Ejector, 3-Way Syringe, 1 Empty			
12	2nd Assi (Hands free suction)			
13	Tissue B	ox Holder		
14	Foot Control			
15	24" Monitor Arm (Mount Type Only)	43" Monitor Rotary Arm (Cart Type Only)		
16	Surgical Table (Cart Type Only)		- Optional	
17	Clean Water System (Distilled Water Bottle)			
18	Spittoon Valve *no drain pipe			



	19	Delivery Table	
	20	Delivery Table (5 Instrument Holders): 3-Way Syringe, 2 High Turbine Tubing	
	21	Micro NX Electric Motor Set (Instrument 3)	Included
	22	Dmetec Scaler - Compact S	
	23	Secondary Tray	
	24	Mouse Table	
=/(25	Micro NX Electric Motor Set (Instrument 2)	Optional

Stool



	26	Functional Doctor Stool	Included
	27	Functional Assistant Stool	Optional

Chair Warranty Coverage

6 Year Warranty:

· Hydraulic system; Motor pump, Solenoid valve.

3 Year Warranty:

· All chairs; Delivery system, Unit, Dental light, Foot control.

2 Year Warranty:

· Doctor stool, Assistant stool, Chair cylinders.

1 Year Warranty:

· All instruments: 3 Way-syringe, Optic tubing, Electric motor set, Scaler set, HVE, SVE.

0 Warranty:

· 3 way-syringe tip, vacuum filter, spittoon filter, damage due to customer's negligence, Seat color changes or tears.

Scope of Warranty

- · This document supersedes any other warranties, agreements, conditions, or other responsibilities, express or implied.
- · If Mega'Gen determines that the product has a defect that is covered under this Warranty, the product will be repaired or replaced with a product that is comparable to the original in performance.
- · If the product is covered under the Mega'Gen Warranty, no part costs will be charged.
- · All product warranty coverage is from date of shipment from Mega'Gen to the Buyer.
- · Products not manufactured by Mega'Gen are covered by the original manufacturer's warranty as specified in the warranty period.
- · All Mega'Gen warranty terms and agreements are non-transferable and apply only to the original device and the original buyer.
- This warranty is void in the event of damage, failure, or defect due to improper maintenance, alteration, or repair by Mega'Gen or anyone other than an authorized dealer.
- The style and color options available for service may vary, and if the product or color is discontinued, a replacement that is most like the existing product may be provided.
- · It is guaranteed to be one year after the part is replaced or larger than the rest of the original warranty period.
- If you have any questions that are not covered by this warranty, please contact Mega'Gen Customer Service.





www.megagenamerica.com

Tel. (844) 288-5425

Email. info@megagenamerica.com